


Medtronic

Designed
to be
durable.

People who know,
think **Evolut™** First.



Designed
to be
durable.

Valve
design impacts
durability.

Durability
impacts
mortality.

Durability starts with design



Built on a proven foundation

With its supra-annular, self-expanding valve frame, Evolut™ TAVR is built on the original CoreValve™ platform which has consistently shown strong EOAs and low gradients over time.

How did we design for durability?

More surface

Taller leaflet mounting allows for a greater distance between the commissure and the edge of the leaflet, distributing stress over a greater distance.


More height

By decoupling the native annular plane where the sealing occurs, from the working portion of the prosthetic leaflets, you can facilitate circularity and maximize leaflet coaptation.

More room

The tall valve keeps the working portion above and unconstrained by the native annulus, allowing for a large effective orifice area.

Supra-annular design benefits



Large EOAs mean less restriction of blood through the valve.

Less restriction leads to low gradients (mean systolic gradient).

Large EOAs have been correlated to less patient-prosthesis mismatch (PPM).

Less PPM and low gradients after aortic valve replacement have been linked to:

- Better survival^{1,2}
- Less heart failure rehospitalization^{2,3}
- Better valve durability^{4,5}

CoreValve™/Evolut™ TAVR platform

Intermediate risk⁶

Average EOA at 5 years (cm²)

Devices used:
83.8% CoreValve
16.2% Evolut™ R

CoreValve/Evolut TAVR platform

Low risk⁷

Average EOA at 2 years (cm²)

Devices used:
3.6% CoreValve
74.1% Evolut R
22.3% Evolut™ PRO

Consistently
strong EOAs

¹ Playford D, et al. *J Am Soc Echocardiogr.* 2020;33:1077-1086.e1.

² Herrmann HC, et al. *J Am Coll Cardiol.* 2018;72:2701-2711.

³ Anand V, et al. *Am J Cardiol.* 2020;125:941-947.

⁴ O'Hair D. Presented at American College of Cardiology 70th Annual Scientific Session & Expo. May 2021.

⁵ Søndergaard L, et al. *J Am Coll Cardiol.* 2019;73:546-553.


⁶ Van Mieghem, et al. 5-Year Clinical and Echocardiographic Outcomes from the Randomized SURTAVI Trial. Presented at TCT 2021.

⁷ Forrest JK, on behalf of the Evolut Low Risk Investigators. The Evolut Low Risk Trial Complete 2-year Follow-up. Presented at EuroPCR 2021.



Supra-annular design benefits

Design elements that produce blood flow stasis and extended blood residence time on the leaflets could increase the risk of thrombosis, resulting in sub-optimal clinical results.¹



Subclinical leaflet thrombosis after TAVR: risk factors, effect on outcome, and treatment options²


RR: Relative risk

SAPT: Single antiplatelet therapy

DAPT: Dual antiplatelet therapy

OAC: Oral anticoagulation

SLT: Subclinical leaflet thrombosis




¹ Midha PA, et al. *Circulation*. 2017;136:1598-1609.

² Bogyi M, et al. *JACC Cardiovasc Interv*. 2021;14:2643-2656.

NOTION¹ 8 years

SVD out to 8 years¹



SAVR	135	113	105	97	84	75	62	54	30
TAVI	139	130	126	115	107	94	80	68	44

The NOTION trial is a multicenter, randomized, head-to-head comparison of CoreValve TAVR versus SAVR followed out to 8 years in lower surgical risk patients ≥ 70 years of age who are eligible for surgery. TAVR had significantly less hemodynamic SVD out to 8 years.

The NOTION 8-year data demonstrates excellent SVD rates in a lower surgical risk patient population. Perhaps most importantly, the data provides a signal of durability for the CoreValve platform versus SAVR.

¹ Søndergaard L. Long-term follow-up of transcatheter and surgical bioprosthetic aortic valves in patients with severe aortic stenosis and lower surgical risk. Presented at PCR Valvese-Course. November 24, 2020.

The CoreValve™ platform was more durable than SAVR at eight years.

> SVD definition

Device used:
100% CoreValve


Valve durability for supra-annular, self-expandable TAV found to be statistically better at five years versus both SAVR and balloon-expandable TAV.

Dr. Attizzani 5-year meta analysis¹

Structural valve deterioration[†]


Only SE performs better than SAVR

Comparison: others versus SAVR (random effects model)



SE performs better than SAVR and BE

Comparison: others versus self-expandable (random effects model)



At five years, supra-annular, self-expandable (SE) valves demonstrated:

- Lowest risk of structural valve deterioration (SVD) compared with balloon-expandable (BE) valves and SAVR.
- Significantly stronger hemodynamics with larger EOAs and lower mean gradients versus BE valves.

Study design

- Meta-analysis
- 10 randomized controlled trials
- 9,388 patients
- Follow-up 1 to 6 years
- Multiple devices[‡]

[†]Based on the longest available follow-up for each of the 10 studies used for this meta-analysis. SVD was defined by the respective authors of each paper.

[‡]CoreValve™, Evolut™ R, Evolut™ PRO, Sapien™*, Sapien 3, Sapien XT, and ACURATE neo™*.

¹ Ueyama H, et al. *Am J Cardiol.* 2021;158:104-111.




CoreValve™
and Evolut™ are
the first and only
TAVR platforms
to demonstrate
a lower SVD
than SAVR.

SVD definition >

Devices used:
88.5% CoreValve
11.5% Evolut™ R

CoreValve and Evolut pooled analysis:

5-year SVD adjusted for competing risk of mortality¹



¹ Reardon et al. 5-Year Incidence, Timing and Predictors of Structural Valve Deterioration of Transcatheter and Surgical Aortic Bioprostheses: Insights from the CoreValve US Pivotal and SURTAVI Trials. Presented at ACC 2022. Updated data on file.













Patients with SVD had a near **two-fold increased risk** for all-cause mortality ($P < 0.001$) and hospitalization for AV disease or worsening heart failure ($P = 0.01$) at five years.

SVD definition 

RCT and Non-RCT cohorts:
97% CoreValve
3% Evolut R

CoreValve™ and Evolut™ pooled analysis:

Worsened clinical outcomes in patients who develop SVD¹

		HR (95% CI)	P value
Pooled surgery RCT and all CoreValve/Evolut (N = 4,762)			
All-cause mortality		2.03 (1.46, 2.82)	< 0.001
Cardiovascular mortality		1.86 (1.20, 2.90)	0.006
Aortic valve-related hospitalization		2.17 (1.23, 3.84)	0.008
Composite [†]		2.02 (1.42, 2.88)	< 0.001
Surgery RCT (N = 971)			
All-cause mortality		2.45 (1.40, 4.30)	0.002
Cardiovascular mortality		2.37 (1.10, 5.08)	0.003
Aortic valve-related hospitalization		2.20 (0.81, 5.98)	0.120
Composite [†]		2.73 (1.53, 4.88)	< 0.001
All CoreValve/Evolut TAVR (N = 3,791)			
All-cause mortality		2.34 (1.55, 3.53)	< 0.001
Cardiovascular mortality		2.17 (1.26, 3.76)	0.006
Aortic valve-related hospitalization		2.45 (1.22, 4.93)	0.010
Composite [†]		2.03 (1.29, 3.19)	0.002


0.10 1.00 10.00
Lower risk with SVD ← → Higher risk with SVD

[†]All-cause mortality or aortic valve-related hospitalization.

¹ Reardon, et al. 5-Year Incidence, Timing and Predictors of Structural Valve Deterioration of Transcatheter and Surgical Aortic Bioprostheses: Insights from the CoreValve US Pivotal and SURTAVI Trials. Presented at ACC 2022. Updated data on file.

Designed to be durable.


1



Established failure rates

NOTION suggested the CoreValve™ platform fails at half the rate of surgery in low-risk patients.


2



Established difference among platforms at five years

Dr. Attizzani established that self-expandable valves demonstrated the lowest risk of SVD compared to balloon-expandable valves and SAVR.

3




Consequence of failure

Dr. Reardon's pooled analysis shows the same statistical trend in durability of SEV over SAVR, as well as the consequence of developing SVD.

Longevity after surgical aortic valve replacement.

Stratification by age and surgical risk groups

Lifetime management of patients undergoing AVR¹



¹ Martinsson A, et al. *J Am Coll Cardiol.* 2021;78:2147-2157.


Supporting data



CoreValve™
TAV remained
hemodynamically
stable at
five years.

CHOICE¹ 5 years

Hemodynamics to 5 years¹



SVD definition >

Device used:
100% CoreValve

For EOAs:

Baseline: $p = 0.71$
Post-TAVR: $p = 0.86$
30 days: $p = 0.13$
1 year: $p = 0.34$
5 years: $p = 0.02$

For gradients:

Baseline: $p = 0.90$
Post-TAVR: $p < 0.001$
30 days: $p < 0.001$
1 year: $p = 0.007$
5 years: $p = 0.001$

In this prospective, randomized study, CoreValve TAV remained hemodynamically stable at 5 years whereas the SAPIEN™ TAV had a 20% decline in EOA and a 40% increase in gradient.

CoreValve also had a statistically significant advantage in terms of freedom from SVD over SAPIEN (0.0% vs. 6.6%; $p = 0.018$).


¹ Abdel-Wahab M, et al. Five-year outcomes after TAVI with balloon-expandable vs. self-expanding valves: Results from the CHOICE randomised clinical trial. Presented at EuroPCR 2019. Paris, France.

Freedom from SVD:
82.4%
 for CoreValve™ TAV
 at seven years.

SVD definition >

Device used:
 100% CoreValve

DEUTSCH¹ 7 years



Retrospective analysis from a single-center registry


This chart clearly demonstrates significantly less SVD for CoreValve than SAPIEN out to 7 years. Freedom from SVD: 82.4% for CoreValve; 63.4% for SAPIEN.

When looking at freedom from SVD, at every time point (1, 3, 5, and 7 years), there was numerically less SVD with CoreValve than with SAPIEN.

¹ Deutsch MA, et al. *EuroIntervention*. 2018;14:41-49.

ITALIAN REGISTRY¹ 8 years

Mean gradient to 8 years¹



Multicenter registry

Together with NOTION, this is the long-term data on the self-expanding, supra-annular CoreValve platform. Data demonstrates very low rates of moderate and severe hemodynamic SVD. The cumulative incidence of moderate and severe SVD at 8 years are 3.0% and 1.6%, respectively.

Additionally, the bioprosthetic valve failure (BVF) was also very low at 2.5% (includes any valve intervention, severe SVD, and any valve-related deaths), signaling durability for the CoreValve platform. The mean gradients remained low through 8 years.

Long-term data on the self-expanding, supra-annular CoreValveTM platform.

> SVD definition

Device used:
100% CoreValve

¹ Testa L, et al. Valve Performance and echocardiographic data throughout 8 years follow up after TAVR. Presented at EuroPCR 2019. Paris, France.



Indications

The Medtronic CoreValve™ Evolut™ R, CoreValve™ Evolut™ PRO, and Evolut™ PRO+ systems are indicated for relief of aortic stenosis in patients with symptomatic heart disease due to severe native calcific aortic stenosis who are judged by a heart team, including a cardiac surgeon, to be appropriate for the transcatheter heart valve replacement therapy.

The Medtronic CoreValve Evolut R, CoreValve Evolut PRO, and Evolut PRO+ systems are indicated for use in patients with symptomatic heart disease due to failure (stenosed, insufficient, or combined) of a surgical bioprosthetic aortic valve who are judged by a heart team, including a cardiac surgeon, to be at high or greater risk for open surgical therapy (e.g., STS predicted risk of operative mortality score ≥ 8% or at a ≥ 15% risk of mortality at 30 days).

Contraindications

The CoreValve Evolut R, CoreValve Evolut PRO, and Evolut PRO+ systems are contraindicated in patients who cannot tolerate Nitinol (titanium or nickel), an anticoagulation/antiplatelet regimen, or who have active bacterial endocarditis or other active infections.

Warnings

General Implantation of the CoreValve Evolut R, PRO, and PRO+ systems should be performed only by physicians who have received Medtronic CoreValve Evolut R, PRO, or PRO+ training. This procedure should only be performed where emergency aortic valve surgery can be performed promptly. Mechanical failure of the delivery catheter system and/or accessories may result in patient complications. *Transcatheter aortic valve (bioprosthesis)* Accelerated deterioration due to calcific degeneration of the bioprostheses may occur in: children, adolescents, or young adults; patients with altered calcium metabolism (e.g., chronic renal failure or hyperthyroidism).

Precautions

General Clinical long-term durability has not been established for the bioprosthesis. Evaluate bioprosthesis performance as needed during patient follow-up. The safety and effectiveness of the CoreValve Evolut R, PRO, and PRO+ systems have not been evaluated in the pediatric population. The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in the following patient populations: Patients who do not meet the criteria for symptomatic severe native aortic stenosis as defined: (1) symptomatic severe high-gradient aortic stenosis – aortic valve area ≤ 1.0 cm² or aortic valve area index ≤ 0.6 cm²/m², a mean aortic valve gradient ≥ 40 mm Hg, or a peak aortic-jet velocity ≥ 4.0 m/s; (2) symptomatic severe low-flow, low-gradient aortic stenosis – aortic valve area ≤ 1.0 cm² or aortic valve area index ≤ 0.6 cm²/m², a mean aortic valve gradient < 40 mm Hg, and a peak aortic-jet velocity < 4.0 m/s; with untreated, clinically significant coronary artery disease requiring revascularization; with a preexisting prosthetic heart valve with a rigid support structure in either the mitral or pulmonic position if either the preexisting prosthetic heart valve could affect the implantation or function of the bioprosthesis or the implantation of the bioprosthesis could affect the function of the preexisting prosthetic heart valve; patients with liver failure (Child-Pugh Class C); with cardiogenic shock manifested by low cardiac output, vasopressor dependence, or mechanical hemodynamic support; patients who are pregnant or breastfeeding. The safety and effectiveness of a CoreValve Evolut R, Evolut PRO, or Evolut PRO+ bioprosthesis implanted within a failed preexisting transcatheter bioprosthesis has not been demonstrated. Implanting a CoreValve Evolut R, Evolut PRO, or Evolut PRO+ bioprosthesis in a degenerated surgical bioprosthetic valve (transcatheter aortic valve in surgical aortic valve [TAV-in-SAV]) should be avoided in the following conditions: The degenerated surgical bioprosthetic valve presents with: a significant concomitant paravalvular leak (between the prosthesis and the native annulus), is not securely fixed in the native annulus, or is not structurally intact (e.g., wire form frame fracture); partially detached leaflet that in the aortic position may obstruct a coronary ostium; stent frame with a manufacturer-labeled inner diameter < 17 mm. The safety and effectiveness of the bioprosthesis for aortic valve replacement have not been evaluated in patient populations presenting with the following: Blood dyscrasias as defined as leukopenia (WBC < 1,000 cells/mm³), thrombocytopenia (platelet count < 50,000 cells/mm³), history of bleeding diathesis or coagulopathy, or hypercoagulable states; congenital unicuspid valve; mixed aortic valve disease (aortic stenosis and aortic regurgitation with predominant aortic regurgitation [3-4+]); moderate to severe (3-4+) or severe (4+) mitral or severe (4+) tricuspid regurgitation; hypertrophic obstructive cardiomyopathy; new or untreated echocardiographic evidence of intracardiac mass, thrombus, or vegetation; native aortic annulus size < 18 mm or > 30 mm for Evolut R/Evolut PRO+ and < 18 mm or > 26 mm for CoreValve Evolut PRO per the baseline diagnostic imaging or surgical bioprosthetic aortic annulus size < 17 mm or > 30 mm for CoreValve Evolut R/Evolut PRO+ and < 17 mm or > 26 mm for Evolut PRO; transarterial access unable to accommodate an 18 Fr sheath or the 14 Fr equivalent EnVeo InLine™ sheath when using Model ENVEOR-US/ENVPRO-14-US/D-EVPROP2329US or transarterial access unable to accommodate a 20 Fr introducer sheath or the 16 Fr equivalent EnVeo InLine sheath when using Model ENVEOR-N-US/ENVPRO-16-US or transarterial access unable to accommodate a 22 Fr introducer sheath or the 18 Fr equivalent Evolut PRO+ InLine sheath when using Model D-EVPROP34US; prohibitive left ventricular outflow tract calcification; sinus of Valsalva anatomy that would prevent adequate coronary perfusion; significant aortopathy requiring ascending aortic replacement; moderate to severe mitral stenosis; severe ventricular dysfunction with left ventricular ejection fraction (LVEF) < 20%; symptomatic carotid or vertebral artery disease; and severe basal septal hypertrophy with an outflow gradient.

Before Use Exposure to glutaraldehyde may cause irritation of the skin, eyes, nose, and throat. Avoid prolonged or repeated exposure to the vapors. Damage may result from forceful handling of the catheter. Prevent kinking of the catheter when removing it from the packaging. The bioprosthesis size must be appropriate to fit the patient’s anatomy. Proper sizing of the devices is the responsibility of the physician. Refer to the Instructions for Use for available sizes. Failure to implant a device within the sizing matrix could lead to adverse effects such as those listed below. Patients must present with transarterial access vessel diameters of ≥ 5 mm when using Model ENVEOR-US/ENVPRO-14-US/D-EVPROP2329US or ≥ 5.5 mm when using Model ENVEOR-N-US/ENVPRO-16-US or ≥ 6 mm when using Model D-EVPROP34US, or patients must present with an ascending aortic (direct aortic) access site ≥ 60 mm from the basal plane for both systems. Implantation of the bioprosthesis should be avoided in patients with aortic root angulation (angle between plane of aortic valve annulus and horizontal plane/vertebrae) of > 30° for right subclavian/axillary access or > 70° for femoral and left subclavian/axillary access. For subclavian access, patients with a patent left internal mammary artery (LIMA) graft must present with access vessel diameters that are either ≥ 5.5 mm when using Models ENVPRO-14-US/ENVEOR-L-US/D-EVPROP2329US or ≥ 6 mm when using Models ENVPRO-16-US and ENVEOR-N-US or ≥ 6.5 mm when using Model D-EVPROP34US. Use caution when using the subclavian/axillary approach in patients with a patent LIMA graft or patent RIMA graft. For direct aortic access, ensure the access site and trajectory are free of patent RIMA or a preexisting patent RIMA graft. For transfemoral access, use caution in patients who present with multiplanar curvature of the aorta, acute angulation of the aortic arch, an ascending aortic aneurysm, or severe calcification in the aorta and/or vasculature. If ≥ 2 of these factors are present, consider an alternative access route to prevent vascular complications. Limited clinical data are available for transcatheter aortic valve replacement in patients with a congenital bicuspid aortic valve who are deemed to be at low surgical risk. Anatomical characteristics should be considered when using the valve in this population. In addition, patient age should be considered as long-term durability of the valve has not been established.

During Use After the procedure, administer appropriate antibiotic prophylaxis as needed for patients at risk for prosthetic valve infection and endocarditis. After the procedure, administer anticoagulation and/or antiplatelet therapy per physician/clinical judgment. Excessive contrast media may cause renal failure. Prior to the procedure, measure the patient’s creatinine level. During the procedure, monitor contrast media usage. Conduct the procedure under fluoroscopy. Fluoroscopic procedures are associated with the risk of radiation damage to the skin, which may be painful, disfiguring, and long-term. The safety and efficacy of a CoreValve Evolut R, Evolut PRO, or Evolut PRO+ bioprosthesis implanted within a transcatheter bioprosthesis have not been demonstrated.

Potential adverse events

Potential risks associated with the implantation of the CoreValve Evolut R, CoreValve Evolut PRO, or Evolut PRO+ transcatheter aortic valve may include, but are not limited to, the following: • death • myocardial infarction, cardiac arrest, cardiogenic shock, or cardiac tamponade • coronary occlusion, obstruction, or vessel spasm (including acute coronary closure) • cardiovascular injury (including rupture, perforation, tissue erosion, or dissection of vessels, ascending aorta trauma, ventricle, myocardium, or valvular structures that may require intervention) • emergent surgical or transcatheter intervention (e.g., coronary artery bypass, heart valve replacement, valve explant, percutaneous coronary intervention [PCI], balloon valvuloplasty) • prosthetic valve dysfunction (regurgitation or stenosis) due to fracture; bending (out-of-round configuration) of the valve frame; underexpansion of the valve frame; calcification; pannus; leaflet wear, tear, prolapse, or retraction; poor valve coaptation; suture breaks or disruption; leaks; mal-sizing (prosthesis-patient mismatch); malposition (either too high or too low)/malplacement • prosthetic valve migration/embolization • prosthetic valve endocarditis • prosthetic valve thrombosis • delivery catheter system malfunction resulting in the need for additional recrossing of the aortic valve and prolonged procedural time • delivery catheter system component migration/embolization • stroke (ischemic or hemorrhagic), transient ischemic attack (TIA), or other neurological deficits • individual organ (e.g., cardiac, respiratory, renal [including acute kidney failure]) or multi-organ insufficiency or failure • major or minor bleeding that may require transfusion or intervention (including life-threatening or disabling bleeding) • vascular access-related complications (e.g., dissection, perforation, pain, bleeding, hematoma, pseudoaneurysm, irreversible nerve injury, compartment syndrome, arteriovenous fistula, or stenosis) • mitral valve regurgitation or injury • conduction system disturbances (e.g., atrioventricular node block, left bundle-branch block, asystole), which may require a permanent pacemaker • infection (including septicemia) • hypotension or hypertension • hemolysis • peripheral ischemia • General surgical risks applicable to transcatheter aortic valve implantation: • bowel ischemia • abnormal lab values (including electrolyte imbalance) • allergic reaction to antiplatelet agents, contrast medium, or anesthesia • exposure to radiation through fluoroscopy and angiography • permanent disability.

Please reference the CoreValve Evolut R, CoreValve Evolut PRO, and Evolut PRO+ Instructions for Use for more information regarding indications, warnings, precautions, and potential adverse events.

Caution: Federal Law (USA) restricts these devices to the sale by or on the order of a physician.

The commercial name of the Evolut™ R device is Medtronic CoreValve™ Evolut™ R System, the commercial name of the Evolut™ PRO device is Medtronic CoreValve™ Evolut™ PRO System, and the commercial name of the Evolut™ PRO+ device is Medtronic Evolut™ PRO+ System.

Medtronic

710 Medtronic Parkway
Minneapolis, MN 55432-5604
USA

Toll-free: 800.328.2518
Tel: +1.763.514.4000

LifeLine
CardioVascular Technical Support
Toll-free: 877.526.7890
Tel: +1.763.526.7890
rs.structuralheart@medtronic.com

medtronic.com

©2022 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic. ™Third-party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.


UC202300043 EN
06/2022



Designed
to be
durable.

Medtronic

The hemodynamic performance of the supra-annular TAVR platform



8-year data from the NOTION trial¹

In a recent update from the NOTION trial on patients at lower surgical risk over the age of 70, we continue to see statistically better hemodynamics versus surgery – as well as a promising signal of durability.


Findings from 8-year follow-up of the NOTION trial¹

- CoreValve™ hemodynamics were statistically better versus surgery out to 8 years.
- CoreValve structural valve deterioration (SVD) was statistically better than surgery out to 8 years.

Structural valve deterioration (SVD)²

- Moderate or greater hemodynamic SVD
- Mean gradient ≥ 20 mm Hg or
- Mean gradient ≥ 10 mm Hg change from baseline or
- Moderate/severe intra-prosthetic aortic regurgitation (AR) (new or worsening from baseline)

Half the rate of SVD versus surgery SVD out to 8 years



SAVR
data yet.

the only platforms
AVR at five years.¹¹
our patients.

4.38%
2.20%
Medtronic TAVR
platforms demonstrated
significantly lower rates
of structural valve
deterioration (SVD)¹
vs. SAVR at five years.

f failure
pooled analysis
statistical trend
EV over SAVR,
consequence of

Designed
to be
durable.

Evolut™ first

Medtronic

EVIDENCE UPDATE


Valve durability

Valve durability for supra-annular, self-expandable TAV found to be statistically better at five years versus both SAVR and balloon-expandable TAV

Structural valve deterioration[†]


Only SE performs better than SAVR

Comparison: others versus SAVR (random effects model)




SE performs better than SAVR and BE

Comparison: others versus self-expandable (random effects model)



Both SE and SAVR perform better than BE

Comparison: others versus balloon-expandable (random effects model)



Study design

- Meta-analysis
- 10 randomized controlled trials
- 9,388 patients
- Follow-up 1 to 6 years
- Multiple devices[‡]

Key observations from the five-year meta-analysis:

At five years, supra-annular, self-expandable (SE) valves demonstrated:

- Lowest risk of structural valve deterioration (SVD) compared with balloon-expandable (BE) valves and SAVR.
- Significantly stronger hemodynamics with larger EOAs and lower mean gradients versus BE valves.

Authors noted that additional studies including newer generations of valves are warranted to address known THV-specific risks, such as AR and reintervention.

SVD was less frequent in SE-THV compared with BE-THV and SAVR (HR 0.14, 95% CI 0.07 to 0.27; HR 0.34, 95% CI 0.24 to 0.47, respectively).

[†]Based on the longest available follow-up for each of the 10 studies used for this meta-analysis. SVD was defined by the respective authors of each paper.

[‡]CoreValve™, Evolut™ R, Evolut™ PRO, Sapien™, Sapien 3, Sapien XT, and ACURATE neo™.



SAVR
data yet.

the only platforms
AVR at five years.¹¹
ur patients.

4.38%
2.20%
Medtronic TAVR
platforms demonstrated
significantly lower rates
of structural valve
deterioration (SVD)[†]
vs. SAVR at five years.

f failure
ooled analysis
statistical trend
EV over SAVR,
nsequence of



Designed
to be
durable.


Evolut™ first

Medtronic

The best TAVR vs. SAVR durability data yet.

CoreValve™ and Evolut™ TAVR systems are the only platforms to demonstrate a durability benefit over SAVR at five years.^{†1}

Valve durability is important to you and your patients.
It's important to us too.



4.38%

2.20%

Medtronic TAVR
platforms demonstrated
significantly lower rates
of structural valve
deterioration (SVD)[‡]
vs. SAVR at five years.


TAVR risks may include, but are not limited to, death, stroke, damage to the arteries, bleeding, and need for permanent pacemaker.

^{†1}In pooled analysis of intermediate and high-risk patients. Devices used: CoreValve 88%/Evolut R 12%.

[‡]Structural valve deterioration (SVD) was defined as an increase in mean gradient ≥ 10 mm Hg over five years with a mean gradient ≥ 20 mm Hg at last echo OR new onset/increase of central AR of \geq moderate in severity.

NOTION¹ 8 years

SVD out to 8 years¹



SVD definition¹

- Moderate or greater hemodynamic SVD
- Mean gradient ≥ 20 mm Hg OR
- Mean gradient ≥ 10 mm Hg change from baseline OR
- Moderate/severe intra-prosthetic aortic regurgitation (AR) (new or worsening from baseline)

¹ Capodanno D, et al. *Eur Heart J*. 2017;38:3382-3390.

SAVR	135	113	105	97	84	75	62	54	30
TAVI	139	130	126	115	107	94	80	68	44

The NOTION trial is a multicenter, randomized, head-to-head comparison of CoreValve TAVR versus SAVR followed out to 8 years in lower surgical risk patients ≥ 70 years of age who are eligible for surgery. TAVR had significantly less hemodynamic SVD out to 8 years.

The NOTION 8-year data demonstrates excellent SVD rates in a lower surgical risk patient population. Perhaps most importantly, the data provides a signal of durability for the CoreValve platform versus SAVR.

¹ Søndergaard L. Long-term follow-up of transcatheter and surgical bioprosthetic aortic valves in patients with severe aortic stenosis and lower surgical risk. Presented at PCR Valvese-Course; November 24, 2020.

CoreValveTM
arm was
durable
AVR at
years.

> SVD definition

Device used:
100% CoreValve

CoreValve[™]
and Evolut[™]
the first an
TAVR platfo
to demons
a lower SV
than SAVR.


CoreValve and Evolut pooled analysis:

5-year SVD adjusted for competing risk of mortality¹

SVD definition¹

SVD was defined as \geq moderate hemodynamic valve deterioration (HVD): Increase in mean gradient \geq 10 mm Hg from discharge/30-day echo to last available echo AND mean gradient \geq 20 mm Hg at last available echo OR new onset/increase of intra-prosthetic aortic regurgitation (AR) \geq moderate.

¹ Adapted from VARC-3 Writing Committee, et al. *Eur Heart J.* 2021;42:1825-1857.



SVD definition >

Devices used:
88.5% CoreValve
11.5% Evolut R

¹ Reardon et al. 5-Year Incidence, Timing and Predictors of Structural Valve Deterioration of Transcatheter and Surgical Aortic Bioprostheses: Insights from the CoreValve US Pivotal and SURTAVI Trials. Presented at ACC 2022. Updated data on file.

CoreValve™ and Evolut™ pooled analysis:

Worsened clinical outcomes in patients who develop SVD¹

Patients with SVD had a near two-fold increased risk of all-cause mortality (P < 0.001), hospitalization for heart failure (P < 0.001) or worse heart failure (P < 0.001) at five years.

SVD definition¹

SVD was defined as ≥ moderate hemodynamic valve deterioration (HVD): Increase in mean gradient ≥ 10 mm Hg from discharge/30-day echo to last available echo AND mean gradient ≥ 20 mm Hg at last available echo OR new onset/increase of intra-prosthetic aortic regurgitation (AR) ≥ moderate.

¹ Adapted from VARC-3 Writing Committee, et al. *Eur Heart J.* 2021;42:1825-1857.




	HR (95% CI)	P value
2)		
	2.03 (1.46, 2.82)	< 0.001
	1.86 (1.20, 2.90)	0.006
	2.17 (1.23, 3.84)	0.008
	2.02 (1.42, 2.88)	< 0.001
	2.45 (1.40, 4.30)	0.002
	2.37 (1.10, 5.08)	0.003
	2.20 (0.81, 5.98)	0.120
	2.73 (1.53, 4.88)	< 0.001
	2.34 (1.55, 3.53)	< 0.001
	2.17 (1.26, 3.76)	0.006
	2.45 (1.22, 4.93)	0.010
	2.03 (1.29, 3.19)	0.002

0.10 1.00 10.00
Lower risk with SVD ← → Higher risk with SVD

SVD definition >

RCT and Non-RCT cohorts:
97% CoreValve
3% Evolut R

[†]All-cause mortality or aortic valve-related hospitalization.
¹ Reardon, et al. 5-Year Incidence, Timing and Predictors of Structural Valve Deterioration of Transcatheter and Surgical Aortic Bioprostheses: Insights from the CoreValve US Pivotal and SURTAVI Trials. Presented at ACC 2022. Updated data on file.



CoreValve[™]
TAV remain
hemodyna
stable at
five years.


CHOICE¹ 5 years

Hemodynamics to 5 years¹

SVD definition¹

- Moderate or greater hemodynamic SVD
- Mean gradient ≥ 20 mm Hg OR
- Mean gradient ≥ 10 mm Hg change from baseline OR
- Moderate/severe intra-prosthetic aortic regurgitation (AR) (new or worsening from baseline)

¹ Capodanno D, et al. *Eur Heart J*. 2017;38:3382-3390.



SVD definition >

Device used:
100% CoreValve

For EOAs:

Baseline: $p = 0.71$
Post-TAVR: $p = 0.86$
30 days: $p = 0.13$
1 year: $p = 0.34$
5 years: $p = 0.02$

For gradients:

Baseline: $p = 0.90$
Post-TAVR: $p < 0.001$
30 days: $p < 0.001$
1 year: $p = 0.007$
5 years: $p = 0.001$

In this prospective, randomized study, CoreValve TAV remained hemodynamically stable at 5 years whereas the SAPIEN[™] TAV had a 20% decline in EOA and a 40% increase in gradient.

CoreValve also had a statistically significant advantage in terms of freedom from SVD over SAPIEN (0.0% vs. 6.6%; $p = 0.018$).

¹ Abdel-Wahab M, et al. Five-year outcomes after TAVI with balloon-expandable vs. self-expanding valves: Results from the CHOICE randomised clinical trial. Presented at EuroPCR 2019; Paris, France.


DEUTSCH¹ 7 years

Freedom from SVD¹

SVD definition¹

- Moderate or greater hemodynamic SVD
- Mean gradient ≥ 20 mm Hg OR
- Mean gradient ≥ 10 mm Hg change from baseline OR
- Moderate/severe intra-prosthetic aortic regurgitation (AR) (new or worsening from baseline)

¹ Capodanno D, et al. *Eur Heart J*. 2017;38:3382-3390.



Patients at risk

Time (years)	0	1	3	5	7	9
CoreValve	214	147	121	106	93	75
SAPIEN	86	57	52	43	36	23

Retrospective analysis from a single-center registry

This chart clearly demonstrates significantly less SVD for CoreValve than SAPIEN out to 7 years. Freedom from SVD: 82.4% for CoreValve; 63.4% for SAPIEN.

When looking at freedom from SVD, at every time point (1, 3, 5, and 7 years), there was numerically less SVD with CoreValve than with SAPIEN.

¹ Deutsch MA, et al. *EuroIntervention*. 2018;14:41-49.


Freedom from
82.4%
for CoreValveTM
at seven years

SVD definition >

Device used:
100% CoreValve

ITALIAN REGISTRY¹ 8 years

Mean gradient to 8 years¹



SVD definition¹

- Moderate or greater hemodynamic SVD
- Mean gradient ≥ 20 mm Hg OR
- Mean gradient ≥ 10 mm Hg change from baseline OR
- Moderate/severe intra-prosthetic aortic regurgitation (AR) (new or worsening from baseline)

¹ Capodanno D, et al. *Eur Heart J*. 2017;38:3382-3390.

Multicenter registry

Together with NOTION, this is the longest-term data on the self-expanding, supra-annular CoreValve platform. Data demonstrates very low rates of moderate and severe hemodynamic SVD. The cumulative incidence of moderate and severe SVD at 8 years are 3.0% and 1.6%, respectively.

Additionally, the bioprosthetic valve failure (BVF) was also very low at 2.5% (includes any valve intervention, severe SVD, and any valve-related deaths), signaling durability for the CoreValve platform. The mean gradients remained low through 8 years.

> SVD definition

Device used:
100% CoreValve

¹ Testa L, et al. Valve Performance and echocardiographic data throughout 8 years follow up after TAVR. Presented at EuroPCR 2019. Paris, France.