


ASAHI SASUKE: the next generation of double lumen catheter



Ordering Information

Product	Catalogue No.	Outer Diameter			Inner Diameter		Usable Length	Recommended Guide Wire	Hydrophilic Coating Length
		Tip	Distal Shaft	Proximal Shaft	Tip	Shaft			
ASAHI SASUKE	SA145-33N	0.50mm (1.5F)	0.84mm/1.08mm (2.5F/3.3F)	1.05mm (3.2F)	0.40mm (0.016")	0.43mm (0.017")	145cm	0.36mm (0.014")	38cm



© 2018 ASAHI INTECC CO., LTD.  
"ASAHI" and "SASUKE" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Manufactured by:




ASAHI INTECC CO.,LTD.  
3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan  
Contact phone number: +81-561-86-9101  
<http://www.asahi-intecc.com>

Distributed by:








# Advantages of ASAHI SASUKE




# ASAHI SASUKE: the next generation of double lumen catheter




6.5mm tip to OTW port - short tip negates the need for deep engagement inside branches and CTOs



Tapered soft tip - provides high trackability




Visible exit - know exactly where the wire is exiting



Oval design - more stability against vessel wall  
a) 1.08mm (3.3Fr) x  
b) 0.84mm (2.5Fr)



Double stainless steel core - kink resistance and pushability



Lubricity - L<sup>3</sup> Coating (long lasting lubricity) maintains trackability